



SURGEON-FOCUSED  
SOLUTIONS



NEEDLE HOLDER PORTFOLIO

## ***An Open Letter to Our Valued Customers***

First of all, we thank you for sharing your compliments, complaints and criticisms regarding our products and services. We rely on this constant dialogue with you, our valued customer, to fuel the ongoing improvement of our products, services and processes.

This process, known as continuous improvement, has infiltrated every department at Marina Medical. In fact, there is an ongoing project that we would like to share with you. The purpose of this project is to design, develop and deliver the most surgeon-focused portfolio of products.

Our strategy for this project is simple - deconstruct and reconstruct our portfolio of products by surgical procedure and work with subject matter experts to identify holes and opportunities for improvement within our portfolio. You will gradually be introduced to our redesigned portfolio of products through a series of procedure-specific catalogs. In this catalog, we are proud to unveil our needle holders.

The project of redesigning our needle holder portfolio consisted of multiple trips to Germany to visit the manufacturing facilities, countless hours of dialogue with our engineers and medical advisors and most importantly... a razor-sharp focus on designing, developing and delivering the most surgeon focused portfolio of needle holders.

We invite you to try our new needle holders and please stay tuned for more surgeon-focused solutions from Marina Medical.

Thank you,

***The Marina Medical Team***

# SURGEON-FOCUSED SOLUTIONS

German tungsten carbide inserts precisely and securely grip the needle.

Bevelled edges prevent the suture from catching.



**DESIGN**



**DEVELOP**



**DELIVER**

We focused on our surgeon's requirements while designing, developing and delivering our brand new portfolio of needle holders.

**To design** our portfolio, we identified which needle holders to offer and opportunities for improving those patterns.


**To develop** our portfolio, we started with the highest quality materials and most capable manufacturing processes to ensure a consistent output.

Finally, **to deliver** our portfolio, we committed to always keep our needle holders in inventory and to always make ourselves available to our customers.

These are just a few of the ways that we design, develop and deliver **surgeon-focused solutions**.


Gold-plated handles identify tungsten carbide inserts.

Spring resistant ratchet feels secure in your hand.



  
**0.5mm**

Designed for suture 3/0 and larger.




  
**0.4mm**

Designed for suture 4/0 to 6/0.



  
**0.2mm**

Designed for suture 7/0 to 10/0.



Designed for the smallest suture sizes.


# NEEDLE HOLDERS

## BABY WEBSTER • WEBSTER




- 420-040** Baby Webster Needle Holder, Serrated, TC, 12cm
- 420-042** Baby Webster Needle Holder, Smooth, TC, 12cm
- 420-060** Webster Needle Holder, Serrated, TC, 13cm
- 420-070** Webster Needle Holder, Smooth, TC, 13cm

### Baby Webster




### Webster



## CRILE-WOOD




- 420-100** Crile-Wood Needle Holder, Serrated, TC, 15cm
- 420-100S** Crile-Wood Needle Holder, Smooth, TC, 15cm
- 420-118** Crile-Wood Needle Holder, Serrated, TC, 18cm
- 420-120** Crile-Wood Needle Holder, Serrated, TC, 20cm
- 420-123** Crile-Wood Needle Holder, Serrated, TC, 23cm





**420-075** Halsey Needle Holder, Serrated, TC, 13cm  
**420-080** Halsey Needle Holder, Smooth, TC, 13cm




# NEEDLE HOLDERS

## OLSEN-HEGAR




- 420-912** Olsen-Hegar Needle Holder, Serrated, TC, 12cm
- 420-913** Olsen-Hegar Needle Holder, Smooth, TC, 12cm
- 420-914** Olsen-Hegar Needle Holder, Serrated, TC, 14cm
- 420-917** Olsen-Hegar Needle Holder, Serrated, TC, 17cm
- 420-918** Olsen-Hegar Needle Holder, Serrated, TC, 19cm




## BAUMGARTNER




- 420-090** Baumgartner Needle Holder, Serrated, TC, 14cm






**420-286**  
**420-288**

Vascular Needle Holder, Serrated, TC, 15cm  
Vascular Needle Holder, Serrated, TC, 18cm




**RYDER**



**420-452**  
**420-454**  
**420-456**  
**420-465**  
**420-466**

Ryder Needle Holder, Serrated, TC, 13cm  
Ryder Needle Holder, Serrated, TC, 15cm  
Ryder Needle Holder, Serrated, TC, 18cm  
Ryder Needle Holder, Serrated, TC, 20cm  
Ryder Needle Holder, Serrated, TC, 23cm




# NEEDLE HOLDERS

## GILLIES




**420-213** Gillies Needle Holder With Scissors, Serrated, TC,



## NEIVERT



**420-216** Neivert Needle Holder, Serrated, TC, 13cm



955 SHOTGUN ROAD  
SUNRISE, FL 33326 USA

TEL: 954.924.4418

TOLL-FREE: 800.697.1119

FAX: 954.924.4419

TOLL-FREE FAX: 800.748.2089

[WWW.MARINAMEDICAL.COM](http://WWW.MARINAMEDICAL.COM)

